



Estd. 1964



YSI-102

Automatic
Tissue Processor



www.yorco.co.in

Latest Improved Version Model

The Latest Model incorporates an Automatic Safety Device and locking system besides other features as in UNI-101. The Safety Device prevents the transfer of Tissue from 10th stage to 11th stage or from 11th to 12th stage, in case the wax has not melted protecting the unit and the Tissues.

The auto locking system checks the transfer of tissues from 12th to 1st stage to safeguard against the accidental reprocessing of the same and locks the machine in the last stage. This lock can be released by pressing of a push switch to prepare the unit for fresh use.

- Two numbers of unique Pin up & Pin down type Imported timers for programming delay as well as the main cycle with full flexibility for the user to program as per customer's requirement without any external tool or tackle or any kind of timing disc.
- Main-process timer for controlling the timing of the tissue basket in each position.
- Delay timer-24 Hrs. Delay timer is provided for weekend and routine delay starting adjustable as per user requirements.

Suitable to work on 220V,
single phase, 50 Hz, AC supply.
Capacity : 1 Liter



Manufactured By :
YORK SCIENTIFIC INDUSTRIES PVT. LTD.


8/6E, Site-IV, Sahibabad Industrial Area, Sahibabad Distt: Ghaziabad (U.P.)201010, Ph: +91-120-4741800, +91-9810083057
email: ysipl@yorco.co.in, web: www.yorco.co.in

Optional / Spare Accessories

- Automatic wax bath Special Type.
- Beaker without li.
- Basket Rotor which moves up-down/ rotates the Tissue basket slowly.
- Tissue Capsule Basket S.S. Even perforations.
- Tissue Capsule S.S. 28 × 28 × 9mm with sliding cover.
- Tissue Capsule S.S. 20×20×10mm with sliding cover.
- Tissue Capsule S.S. Circular, 36mm dia. with auto press lid.
- Divider for Tissue capsule 'S' shape.
- Divider for Tissue capsule 'V' shape.

Supplied Complete With

- | | |
|---------------------------------------|--------------|
| ● S.S. Basket Rotor |1 No. |
| ● S.S. Tissue Basket |1 No. |
| ● S.S. Tissue Capsule |24 Nos. |
| ● Glass Beaker of 1 lt. Capacity each |10 Nos. |
| ● Beaker Cover |11 Nos. |
| ● S.S. Wax Bath |2 Nos. |



Manufactured By :
YORK SCIENTIFIC INDUSTRIES PVT. LTD.

8/6E, Site-IV, Sahibabad Industrial Aea, Sahibabad Distt: Ghaziabad (U.P.)201010, Ph: +91-120-4741800, +91-9810083057
email: ysipl@yorco.co.in, web: www.yorco.co.in